

### graphing the black body and visual response spectra

$T_{\text{color}} \equiv 2856\text{-K}$  color temperature of source used for CCD measurement


$\text{vis} := 0..100$   $E\Lambda_{\text{max}} := 700\cdot\text{nano}\cdot\text{m}$   $E\Lambda_{\text{min}} := 400\cdot\text{nano}\cdot\text{m}$

$$\lambda_{\text{vis}} := E\Lambda_{\text{min}} + (E\Lambda_{\text{max}} - E\Lambda_{\text{min}}) \cdot \frac{\text{vis}}{100}$$

$BB_{\text{vis}} := BB_{\lambda}(\lambda_{\text{vis}}, T_{\text{color}})$  normalized black body

$EYE_{\text{vis}} := VIS(\lambda_{\text{vis}})$  photopic vision

$P_{\text{vis}} := BB_{\text{vis}} \cdot EYE_{\text{vis}}$  product of the above



red : black body  
blue : normalized photopic eye response  
green : eye response to black body spectrum

$$VRI := \int_{E\Lambda_{\text{min}}}^{E\Lambda_{\text{max}}} BB_{\lambda}(\lambda, T_{\text{color}}) \cdot VIS(\lambda) d\lambda \quad VRI = 3.724 \times 10^{-8} \text{ m} \quad \begin{array}{l} \text{VRI} \\ \text{integral of} \\ \text{eye response to black body} \\ \text{spectrum} \end{array}$$

$S_v := 100 \cdot I_x$  illumination on imager : specified in imager data sheet

$$K_E := \frac{S_v}{K_p \cdot VRI} \quad K_E = 3.99 \text{ W} \cdot \text{m}^{-2} \cdot \text{micron}^{-1}$$

$$M_{\text{peak}}(T_{\text{color}}) = 2.448 \times 10^6 \text{ W} \cdot \text{m}^{-2} \cdot \text{micron}^{-1}$$

**graphing the detector response**

$$D\lambda_{\min} := 200 \cdot \text{nano} \cdot \text{m} \quad D\lambda_{\max} := 1100 \cdot \text{nano} \cdot \text{m} \quad \text{wavelength range}$$


$$\text{sub} := 0..100 \quad \text{calculation control subscript}$$

$$\lambda D_{\text{sub}} := D\lambda_{\min} + (D\lambda_{\max} - D\lambda_{\min}) \cdot \frac{\text{sub}}{100} \quad \text{wavelength of subth calculation}$$

$$\text{DET}_{\text{sub}} := \text{DSR}(\lambda D_{\text{sub}}) \quad \text{detector spectral response}$$

$$\text{DBB}_{\text{sub}} := \text{BB}_{\lambda}(\lambda D_{\text{sub}}, T_{\text{color}}) \quad \text{black body spectrum}$$

$$\text{PP}_{\text{sub}} := \text{DBB}_{\text{sub}} \cdot \text{DET}_{\text{sub}} \quad \text{detector spectral response to black body excitation}$$



red : black body  
 blue : detector spectral  
 response  
 green : detector response to  
 black body spectrum

**calculate detector experimental responses**

$$\text{DE} := \int_{D\lambda_{\min}}^{D\lambda_{\max}} \text{BB}_{\lambda}(\lambda, T_{\text{color}}) \cdot \text{DSR}(\lambda) \, d\lambda \quad \text{DE} = 1.358 \times 10^{-7} \text{ m} \quad \text{DE detector effectivity}$$



$$\text{DSI}(\lambda) := K_E \cdot \text{BB}_{\lambda}(\lambda, T_{\text{color}}) \quad \text{detector spectral irradiance}$$

$$\text{DI} := \int_{D\lambda_{\min}}^{D\lambda_{\max}} \text{DSI}(\lambda) \, d\lambda \quad \text{DI} = 1.879 \text{ W} \cdot \text{m}^{-2} \quad \text{detector irradiance}$$

$$\text{PSI}(\lambda) := \frac{\text{DSI}(\lambda)}{\left( \frac{h \cdot c}{\lambda} \right)} \quad \text{function for detector spectral photon irradiance}$$

$$\text{PI} := \int_{D\lambda_{\min}}^{D\lambda_{\max}} \text{PSI}(\lambda) \, d\lambda \quad \text{PI} = 7.961 \times 10^{18} \frac{1}{\text{m}^2 \cdot \text{s}} \quad \text{detector photon irradiance}$$

red : black body  
 blue : normalized photopic eye  
 response  
 green : eye response to black  
 body spectrum



**detector geometry**

$$H_{\text{pixel}} := 0.50 \cdot \text{mm} \quad W_{\text{pixel}} := 25 \cdot \text{micron} \quad A_{\text{pixel}} := H_{\text{pixel}} \cdot W_{\text{pixel}} \quad A_{\text{pixel}} = 1.25 \times 10^{-4} \text{ cm}^2$$

$$PI \cdot A_{\text{pixel}} = 9.951 \times 10^{10} \text{ s}^{-1} \quad \text{photon flux on a pixel}$$

$$\text{EPSI}(\lambda) := \text{PSI}(\lambda) \cdot \text{DSR}(\lambda) \quad \text{effective detector spectral photon irradiance}$$

$$\text{EPI} := \int_{D\lambda_{\text{min}}}^{D\lambda_{\text{max}}} A_{\text{pixel}} \text{EPSI}(\lambda) d\lambda \quad \begin{array}{l} \text{effective detector photon flux} \\ \text{these are the photons that are} \\ \text{detected} \end{array} \quad \text{EPI} = 2.311 \times 10^{10} \text{ s}^{-1}$$

$$E_{\text{sat}} := 0.145 \cdot \text{lx} \cdot \text{s} \quad \text{saturation exposure}$$

$$t_{\text{sat}} := \frac{E_{\text{sat}}}{S_v} \quad \text{corresponding exposure time} \quad t_{\text{sat}} = 1.45 \times 10^{-3} \text{ s}$$

$$\text{SPE} := \text{EPI} \cdot t_{\text{sat}} \quad \text{saturation photon exposure} \quad \text{SPE} = 3.351 \times 10^7$$

$$\text{SEC} := \text{SPE} \cdot q \quad \text{saturation charge on pixel} \quad \text{SEC} = 5.369 \times 10^{-12} \text{ C}$$

Professor Jim Palmers Radiometry FAQ page University of Arizona  
<http://www.optics.arizona.edu/Palmer/rpfaq/rpfaq.htm>

**constants and conversions**

$$h \equiv 6.626 \cdot 10^{-34} \text{ J} \cdot \text{s} \quad \text{Planck's constant}$$

$$k \equiv 1.381 \cdot 10^{-23} \text{ J} \cdot \text{K}^{-1} \quad \text{Boltzmann's constant}$$

$$q \equiv 1.602 \cdot 10^{-19} \text{ C} \quad \text{electron charge}$$

$$\text{eV} \equiv 1.602 \cdot 10^{-19} \text{ J} \quad \text{electron volt}$$

$$c \equiv 2.998 \cdot 10^8 \text{ m} \cdot \text{sec}^{-1} \quad \text{speed of light}$$

$$\sigma \equiv 5.670 \cdot 10^{-8} \text{ W} \cdot \text{m}^{-2} \cdot \text{K}^{-4} \quad \text{Stefan-Boltzmann constant}$$

$$K_p \equiv 673 \cdot \text{lm} \cdot \text{W}^{-1} \quad \text{photopic luminous efficacy @ 555 nm}$$

$$\text{milli} \equiv 10^{-3} \quad \text{micro} \equiv 10^{-6} \quad \text{nano} \equiv 10^{-9} \quad \text{pico} \equiv 10^{-12}$$

$$\text{micron} \equiv \text{micro} \cdot \text{m} \quad \text{nm} \equiv \text{m} \cdot 10^{-9}$$

**GLOBAL DEFINITIONS OF FUNCTIONS****black body  
spectral radiant exitance**

$$M_{\lambda}(\lambda, T) \equiv \frac{2 \cdot \pi \cdot h \cdot c^2}{\lambda^5 \cdot \left( e^{\frac{h \cdot c}{\lambda \cdot k \cdot T}} - 1 \right)}$$

$$\text{Typical\_Units} := \text{W} \cdot \text{m}^{-2} \cdot \text{micron}^{-1}$$

**Wein's displacement law**

$$\Lambda_{\text{peak}}(T) \equiv \frac{2898}{T} \cdot \text{micro} \cdot \text{m} \cdot \text{K}$$

This gives the wavelength of the peak of the black body spectrum at temperature T

$$\text{Typical\_Units} := \text{micron}$$

**spectral radiant exitance  
at black body peak**

This gives the spectral radiant exitance of the peak of the black body spectrum at temperature T

$$M_{\text{peak}}(T) \equiv M_{\lambda}(\Lambda_{\text{peak}}(T), T)$$

$$\text{Typical\_Units} := \text{W} \cdot \text{m}^{-2}$$

**normalized black body**

This gives the spectrum, normalized to a peak of 1, of a black body at temperature T

$$\text{BB}_{\lambda}(\lambda, T) \equiv \frac{M_{\lambda}(\lambda, T)}{M_{\text{peak}}(T)}$$

$$\text{Units} := \text{dimensionless}$$

**black body total radiant exitance**

This gives the total radiant power emitted from a black body at temperature T.

$$\text{BBP}(T) \equiv \sigma \cdot T^4$$

$$\text{Typical\_Units} := \text{W} \cdot \text{m}^{-2}$$

**radiant exitance**

This gives the radiant exitance over a passband  
center wavelength =  $\lambda$   
band width =  $\Delta\lambda$   
for a black body spectrum at temperature T

$$M(\lambda, \Delta\lambda, T) \equiv \left( \int_{\lambda - \frac{\Delta\lambda}{2}}^{\lambda + \frac{\Delta\lambda}{2}} M_{\lambda}(\lambda, T) d\lambda \right)$$

$$\text{Typical\_Units} := \text{W} \cdot \text{m}^{-2}$$

## GLOBAL DEFINITION OF PHOTOPIC VISION CURVE

(	400	)			(	.0004	)
	410					.0012	
	420					.0040	
	430					.0116	
	440					.0230	
	450					.0380	
	460					.0600	
	470					.0910	
	480					.1390	
	490					.2080	
	500					.3230	
	510					.503	
	520					.7100	
	530					.8620	
lamda ≡	540	·10 <sup>-9</sup> ·m		photopic ≡		.9540	
	550					.9950	
	560					.9950	
	570					.9520	
	580					.8700	
	590					.7570	
	600					.6310	
	610					.5030	
	620					.3810	
	630					.2650	
	640					.1750	
	650					.1070	
	660					.0610	
	670					.0320	
	680					.0170	)

vision ≡ cspline(lamda, photopic)

VIS(color) ≡ interp(vision, lamda, photopic, color)

VIS(500·nano·m) = 0.323

**GLOBAL DEFINITION OF DETECTOR SPECTRAL RESPONSE CURVE****DETECTOR** Hamamatsu S8377/S8378 series

200	0.24
205	0.245
225	0.29
230	0.32
235	0.29
241	0.25
243	0.24
245	0.25
247	0.295
250	0.31
255	0.295
260	0.28
265	0.27
270	0.28
283	0.465
290	0.48
295	0.465
302	0.44
308	0.43
314	0.44
320	0.645
327	0.67
334	0.645
352	0.53
$\Lambda \equiv 360$	$R \equiv 0.51$
368	0.53
375	0.875
385	0.89
395	0.875
403	0.62
410	0.58
420	0.62

465	0.95
480	1
500	0.95
525	0.84
540	0.81
555	0.84
580	0.9
600	0.93
620	0.9
670	0.6
700	0.5
740	0.4
800	0.3
840	0.2
900	0.1
1000	0.04
1100	0

interpolation functions for the detector spectral response curve

$$\text{Ham} \equiv \text{cspline}(\Lambda, \mathbf{R})$$

$$\text{DSR}(\text{color}) \equiv \text{interp}(\text{Ham}, \Lambda, \mathbf{R}, \text{color}) \quad \text{detector spectral response}$$



